


Note: Part numbers and standards' values may change. Consult LECO for the latest information.



See reverse side for contact information

FUSION SUPPLIES

501-598 Etchable Nibbled Nickel 100 g bottle
501-073 Graphite Powder 100 g bottle
780-892 Crucible Cover 1000/package
502-040 Large Tin Capsule 100/bottle
611-351-180 Lower Electrode Tip (for 764-330 Crucibles with Automation)
611-351-183 Lower Electrode Tip (for 764-330 Crucibles with Automation)
619-896 Lower Electrode Tip (for 619-895 Crucibles)
502-822 Nickel Capsule 100/bottle
763-029 1.5 g Nickel Baskets 100/bottle
502-345 1.5 g Nickel Baskets Ultra High Purity 100/bottle
502-344 1.0 g Nickel Baskets Ultra High Purity 100/bottle
501-059 Small Tin Capsule 100/bottle
761-744 Molecular Sieve Assembly

CATALYSTS AND REAGENTS

501-171-HAZ Anhydron 454 g (1 lb.) bottle
502-189 Copper Sticks 100 g ampule
502-304 Copper Sticks 100 g bottle
502-304-500 Copper Sticks 500 g bottle
502-295 Copper Turnings 200 g bottle
501-081 Glass Wool 454 g (1 lb.) box
501-073 Graphite Powder 100 g bottle
502-174-HAZ Lecosorb 500 g bottle, 20-30 mesh
502-176-HAZ Lecosorb 500 g bottle, 8-20 mesh
761-747-HAZ Schutze™ Reagent

OPERATION ACCESSORIES

766-053 Crucible Tweezers
766-646 Air Filter Coating
775-306 Secondary Filter, 10 micron
775-307 Screen Filter Removal Tool (Thumbscrew)
765-933 Lubricool Cooling Agent
764-078 One gram Weight
783-786-HAZ Oxygen and Moisture Removal Kit
783-785-110-HAZ Replacement Tube for 783-786 Kit
760-138 Sample Tweezers (Straight)
501-241 Vacuum Grease
767-473 Utility Funnel
782-880 Vacuum Cleaner Filter
618-809 Medium Reagent Tube (Glass)

NOTE: Part numbers and standards' values may change. Consult LECO for latest information.

Please Note: Quartz Tubes are typically utilized for higher temperature applications (i.e. catalyst tubes). Glass reagent/filter tubes are designed for lower temperature applications. Glass tubes should not be used in place of quartz tubes.